

# ImageStream® FISHIS® for Blood Cells

## Reagents:

1. Probes: SpectrumGreen Chromosome Enumeration Probes (CEP) for X, Y, or 8 (Abbott Molecular)
2. Carnoy's Fixative (3 parts Methanol : 1 part Glacial Acetic Acid) make fresh, keep cold
3. 2X SSC
4. PBS
5. NP-40

## Protocol:

### 1. Cell fixation:

- a. Enrich mononuclear cells from donor blood by ficoll-hypaque density centrifugation OR use cultured cell line
- b. Pellet PBS-washed cells ( $1-2 \times 10^7$  cells) in 15 ml polystyrene centrifuge tube (centrifuge 300xg for 8 min)
- c. Resuspend in 1.0 ml ice-cold PBS and transfer to microfuge tube. While vortexing, add 0.4 ml ice-cold 100% Carnoys dropwise. Incubate for 10 minutes at room temperature.
- d. Store at -20C for at least 4 hours and up to 3 months

### 2. Hybridization:

- a. Centrifuge 300xg for 8 min. Carefully remove all of the supernatant with pipet.
- b. Resuspend in 1mL of 0.1%NP-40 in 2X SSC buffer.
- c. Aliquot 50 $\mu$ l (~1-2million) cells per hybridization into siliconized microfuge tubes.
- d. Centrifuge 300xg for 8 min. Carefully remove all of the supernatant with pipet.
- e. Per hybridization reaction, prepare a master mix of 28  $\mu$ l hybridization buffer (CEP kit hybridization buffer), 2  $\mu$ l of probe and 10  $\mu$ l of nuclease-free water.
- f. Resuspend cell pellet in 39  $\mu$ l of the master mix and transfer to PCR tube
- g. Thermocycler conditions: 5 min 80°C, 9hr 42°C; optional 4°C storage step.
- h. Add 200 $\mu$ l of 0.1%NP-40 in 2x SSC and transfer to siliconized microfuge tube.
- i. Centrifuge 300xg for 8 min.
- j. Resuspend in 200  $\mu$ l pre-warmed (73°C) 0.3% NP-40/2xSSC.
- k. Incubate at 73°C for 2 minutes, then add 200uL ice-cold PBS
- l. Centrifuge 300xg for 8 min.
- m. Resuspend in 60  $\mu$ L ice-cold PBS and run on ImageStreamX

### 3. Controls:

- a. Monosomy detection: Spike in known amount of CEP-X probed male into CEP-X-probed female.
- b. Trisomy detection: Spike in known amount of CEP-8+Y-probed male into CEP-8-probed male PBL