

# ImageStream® Bacteria Live/Dead Protocol

## Samples: (~1x10<sup>8</sup> bacteria/ml) per sample

Bacteria	Sample 1	Sample 2	Sample 3	Sample 4	Sample 5
<b>Live:Dead</b>	100% to 0%	75% to 25%	50% to 50%	25% to 75%	0% to 100%
<b>Volume</b>	100ulL:0ulD	75ulL:25ulD	50ulL:50ulD	25ulL:75ulD	0ulL:100ulD

## Materials

- 1.) **0.85% NaCl** solution in dH<sub>2</sub>O (Sigma #)
- 2.) **LB-Media** (Sigma #)
- 3.) **70% Isopropanol** (Sigma #)
- 4.) **Propidium Iodide** 20mM Stock (Invitrogen #)
- 5.) **E. Coli** starter Culture (Invitrogen #)

## Bacteria Preparation

- 1.) Add 100ul LB media to sterile erlenmeyer flask using sterile technique.
- 2.) Inoculate with a single colony from E. Coli starter culture.
- 3.) Place on shaker and incubate room temperature overnight.

## Live, Dead Cell Preparation

- 1.) To 2x 1.5ml eppendorf tubes add 1ml of growing culture.
- 2.) Centrifuge 10,000rpm for 5min.
- 3.) Wash 1x in 1ml 0.85% NaCl buffer.
- 4.) Centrifuge 10,000rpm for 5min,
  - a. To one tube resuspend in 0.85%NaCl buffer (this is the live sample)
  - b. To the second: disrupt the pellet by vigorous mixing and resuspend in 100ul of 0.85% NaCl mixing thoroughly to prevent clumps. Next add 1ml of 70% Isopropanol and vortex vigorously (this is the dead sample).
- 5.) Incubate 30min at room temperature mixing every 15min.
- 6.) Centrifuge both samples 10,000rpm for 5min, and resuspend in 1ml 0.85% NaCl.

## Bacteria Staining Protocol

- 1.) Prepare 5x 1.5ml eppendorf tubes with 890ul of 0.85%NaCl.
- 2.) Add 10ul of 20mM (or equivalent) PI.
- 3.) Mix each according to the dilution factors listed in the table below:

Bacteria	Sample 1	Sample 2	Sample 3	Sample 4	Sample 5
<b>Live:Dead</b>	100% to 0%	75% to 25%	50% to 50%	25% to 75%	0% to 100%
<b>0.85% NaCl</b>	890ul	890ul	890ul	890ul	890ul
<b>20mM PI</b>	10ul	10ul	10ul	10ul	10ul
<b>Volume</b>	100ulL:0ulD	75ulL:25ulD	50ulL:50ulD	25ulL:75ulD	0ulL:100ulD
<b>Live:Dead Ratio</b>					

- 4.) Incubate for 15min in the dark at room temperature.
- 5.) Vortex vigorously and aliquot 100ul into a 1.5ml eppendorf tube for running on ImageStream.
- 6.) Dilute 10ul of the 100% Live and 100% Dead samples 1 to 10 with 0.85%NaCl for counting on a hemocytometer.
- 7.) Analyze data in IDEAS and generate a live:dead ratio plot in excel.

# ImageStream® Bacteria Live/Dead Protocol

## Instrument Settings

Load INSPIRE template named "Bacteria PI". To have the fluidics changes load properly a sample must be loaded and running on the system, otherwise an error is thrown and the fluidics portion is entered manually. Key changes to the INSPIRE template are listed below.

- 1.) Turn off un-used channels (collect channels 1, 4, and 6 only for PI).
- 2.) Magnification is 40x.
- 3.) 488 laser power to 100mw (or a level that maximizes PI fluorescence without saturation).
- 4.) 785 laser power to 70mw.
- 5.) 4um core diameter (runs 1ul every ~5min)
- 6.) 10% Beads.
- 7.) 60mm/second core velocity.
- 8.) Turn off auto bead detect.
- 9.) Set squelch to -10.
- 10.) Remove all cell classifiers to collect everything.
- 11.) Collect 50,000 total events to get more than 10,000 bacteria.

## Counting Bacteria With the Hemocytometer

- 1.) To 90ul of 0.85% NaCl add 10ul of bacteria in single cell suspension.
- 2.) Vortex vigoursly
- 3.) Place 10ul under hemocytometer cover slip.
- 4.) Let the sample settle for 20min.
- 5.) Use 40x and count bacteria being aware that some cells may be out of the normal plane of focus and may require some focus adjustment.
- 6.) Count 10 of the smallest squares (these are 1:16<sup>th</sup> of the 1x10<sup>4</sup><sup>th</sup> square).
- 7.) Average all 10 counts.
- 8.) Calculate the final concentration:

**Concentration = Average x dilution x factor x 1x10<sup>4</sup><sup>th</sup> cells per ml**

**Concentration = [Ave x 10 x 16 x 1x10<sup>4</sup><sup>th</sup>] cells/ml**